

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date

Atorvastatin EP Impurity H/Atorvastatin USP Related Compound H/Atorvastatin Lactone

Catalogue no.: VL10002

Batch No.

IUPAC Name: 5-(4-Fluorophenyl)-2-(1-methylethyl)-N,4-diphenyl-1-[2-[(2R,4R)-tetrahydro-4-hydroxy-6-oxo-2H-pyran-2-yl]ethyl]-1H-pyrrole-3-carboxamide



Molecular Formula: C₃₃H₃₃FN₂O₄

Molecular weight: 540.62 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Unknown	-----
2	Purity by HPLC	>95.00 %	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by