

## CERTIFICATE OF ANALYSIS


Reference Substance

Date

**Oxo Atorvastatin Impurity**

Batch No.

IUPAC Name: 4-[6-(4-Fluoro-phenyl)-6-hydroxy-1b-isopropyl-6a-phenyl-1a-phenylcarbamoyl-hexahydro-1,2-dioxo-5a-aza-cyclo-propa[a]inden-3-yl]-3-(R)-hydroxy-butyric acid.


 Molecular Formula: C<sub>33</sub>H<sub>35</sub>FN<sub>2</sub>O<sub>6</sub>

Molecular weight: 574.64g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Light yellow colored	-----
2	Purity by HPLC	>95.00%	Chromatogram attached
3	<sup>1</sup> H NMR	Conform	Report attached
4	MASS	Conform	Report attached

**Storage condition:** Store at 2-8°C, protect from direct sunlight and moisture.

**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:**

**Retest date:**

**Document data reference:**

*Prepared by*

*Reviewed by*

*Approved by*