

## CERTIFICATE OF ANALYSIS


Reference Substance

Date:

**Atorvastatin Related Compound I/ Acetonide Impurity**

Batch No.

IUPAC Name: *Tert*-Butyl 2-((4*R*,6*R*)-6-{2-[2-(4-fluorophenyl)-5-isopropyl-3-phenyl-4-(phenylcarbamoyl)-1*H*-pyrrol-1-yl]ethyl}-2,2-dimethyl-1,3-dioxan-4-yl)acetate.



Molecular Formula: C<sub>40</sub>H<sub>47</sub>FN<sub>2</sub>O<sub>5</sub>

Molecular weight: 654.81

Sr.#	Test	Results	Data ref.
1	Appearance	Off white colored powder	-----
2	Purity by HPLC	>95.0%	Chromatogram attached
3	<sup>1</sup> H NMR	Conform	Report attached
4	MASS	Conform	Report attached

**Storage condition:** Store at 2-8°C, Protect from direct sunlight. Protect from moisture.

**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Date of manufacturing:**

**Retest date:**

**Document data reference:**

*Prepared by*

*Reviewed by*

*Approved by*