

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date

Atorvastatin Diepoxide Calcium Salt

Catalogue no.: VL10014

Batch No.

IUPAC Name: 7-[N-Phenyl-2-[4-(4-fluorophenyl)-6-methylethyl]-3,7-dioxo-5-azatricyclo [4.1.0.0.2,4]-4-[[phenylamino]carboxamide]-(3R,5R)-3,5-dihydroxyheptanoic Acid Calcium Salt



Molecular Formula: C₆₆H₆₈CaF₂N₄O₁₄

Molecular weight: 1219.34 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off-White to Pale Yellow	-----
2	Purity by HPLC	>95.00 %	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by