

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date

4-Hydroxy Atorvastatin Monosodium Salt

Catalogue no.: VL10015

Batch No.

IUPAC Name: Sodium (3R,5R)-7-(2-(4-fluorophenyl)-4-((4-hydroxyphenyl)pyrrol-1-yl)-3,5-dihydroxyheptanoate



Molecular Formula: C₃₃H₃₄FN₂NaO₆

Molecular weight: 596.62 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Unknown	-----
2	Purity by HPLC	>95.00 %	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by