

CERTIFICATE OF ANALYSIS


Reference Substance

Date-

Metformin Related Compound C

Batch No.:

IUPAC Name: N,N-Dimethyl-[1,3,5]triazine-2,4,6-triamine

Molecular Formula: C₅H₁₀N₆

Molecular Weight: 154.17 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White Solid	-----
2	Purity by HPLC	99.96%	chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C and protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:**

Prepared by

Reviewed by

Approved by