

## CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date:

**Metamizole impurity B****Catalogue no.: VL105002**

Batch No.:

IUPAC Name: 4-Aminoantipyrine

Molecular Formula: C<sub>11</sub>H<sub>13</sub>N<sub>3</sub>O

Molecular Weight: 203.25g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White solid	-----
2	Purity by HPLC	>95.00%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

**Storage condition:** Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*