

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date

Leflunomide Impurity A

Catalogue no.:VL1060002

Batch No.

IUPAC Name: 4-(trifluoromethyl) aniline

Molecular Formula: C₇H₆F₃N

Molecular weight: 161.12/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White color sticky solid	-----
2	Purity by HPLC	99.96%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*