


CERTIFICATE OF ANALYSIS

Reference Substance

Date

Moxifloxacin Impurity C

Batch No.

 IUPAC Name: 1-Cyclopropyl-8-ethoxy-6-fluoro-1,4-dihydro-7-[(4*a*S,7*a*S)-octahydro-6*H*-pyrrolo[3,4-*b*]pyridin-6-yl]-4-oxo-3-quinolinecarboxylic Acid

 Molecular Formula: C₂₂H₂₆FN₃O₄

Molecular weight: 415.46g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Yellow solid	-----
2	Purity by HPLC	99.77%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:
Retest date:
Document data reference:
Prepared by
Reviewed by
Approved by