

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Olanzapine Demethyl Impurity

Batch No.

IUPAC Name: 2-Methyl-4-(1-piperazinyl)-10Hthieno[2,3-b][1,5]benzodiazepine



Molecular Formula:C16H18N4S

Molecular weight:298.4 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Light yellow powder	-----
2	Purity by HPLC	99.9%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at Refrigerator, Protect from direct sunlight.**Shipment date:**

This certificate of analysis is valid for one year from date of shipment, provided the substance is stored under suitable conditions.

Retest date:**Document data reference:***Prepared by**Reviewed by**Approved by*