

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Perindopril EP Impurity E

Batch No.:

IUPAC Name: (2S,3aS,7aS)-1-[(2S)-2-[[[(1S)-1-[(1-Methylethoxy)carbonyl]butyl]amino]propanoyl]octahydro-1H-indole-2-carboxylic acid



Molecular Weight: 382.49 g/mol

 Molecular Formula: C₂₀H₃₄N₂O₅

Sr.#	Test	Results	Data ref.
1	Appearance	White solid Powder	-----
2	Purity by HPLC	97.17%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from moisture and direct sunlight.

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:
Retest Date:
Document data reference:
Prepared by
Reviewed by
Approved by