

## CERTIFICATE OF ANALYSIS


Reference Substance

Date

**Phenylephrine Impurity-E**

Batch No.

IUPAC Name: 2-(benzylmethylamino)-1-(3-hydroxyphenyl)ethanone

Molecular Formula: C<sub>16</sub>H<sub>17</sub>NO<sub>2</sub>

Molecular weight: 255.31 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	>95.0%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

**Storage condition:** Store at 2-8°C, Protect from direct sunlight and moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*