


## CERTIFICATE OF ANALYSIS

Reference Substance

Date:

**Benzyl phenylephrine****Batch No.:**

IUPAC Name: (R)-3-(2-Amino-1-hydroxyethyl)phenol

Molecular Formula: C<sub>16</sub>H<sub>19</sub>NO<sub>2</sub>

Molecular Weight: 257.33 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Sticky Mass	-----
2	Purity by HPLC	94.49%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

**Storage condition:** Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*