

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Rosuvastatin Related Compound B/ Rosuvastatin EP Impurity E

Batch No.:

IUPAC Name: (3R,5S,E)-7-(4-(4-fluorophenyl)-2-(2-(4-(4-fluorophenyl)-6-isopropyl-2-(N-methylmethylsulfonamido)pyrimidin-5-yl)-2-hydroxy-N-methylethylsulfonamido)-6-isopropylpyrimidin-5-yl)-3,5-dihydroxyhept-6-enoic acid



Molecular Weight: 831.93 g/mol

 Molecular Formula: C₃₈H₄₅F₂N₆O₉S₂ · 1/2Ca²⁺

Sr.#	Test	Results	Data ref.
1	Appearance	Off White solid	-----
2	Purity by HPLC	95.04%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from moisture and direct sunlight.

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:
Retest Date:
Document data reference:
Prepared by
Reviewed by
Approved by