


CERTIFICATE OF ANALYSIS

Reference Substance
Date:
Rosuvastatin Z-Isomer
Batch No.:
IUPAC Name: (3R,5R,6E)-7-[4-(4-fluorophenyl)-2-(N-methylmethanesulfonamido)-6-(propan-2-yl)pyrimidine-5-yl]-3-dihydroxyhept-6-enoic acid calcium salt

Molecular Formula: C₂₂H₂₈FN₃O₆S

Molecular Weight: 481.54

Sr.#	Test	Results	Data ref.
1	Appearance	Off white Solid	-----
2	Purity by HPLC	99.64%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture. It is moisture sensitive.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by