

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Acyclovir Impurity A / O-Acetyl Acyclovir

Batch No.

IUPAC Name: 2-((2-amino-6-oxo-1H-purin-9(6H)-yl) methoxy) ethyl acetate.

Molecular Formula: C₁₀H₁₃N₅O₄

Molecular weight: 267.24 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off White Solid Powder	-----
2	Purity by HPLC	99.04%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and Moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest date:**

Document data reference:

*Prepared by**Reviewed by**Approved by*