


CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date

Acyclovir EP Impurity F/ N2-Acetyl Acyclovir**Catalogue no.: VL13007**

Batch No.

IUPAC Name: N-[6,9-Dihydro-9-[(2-hydroxyethoxy)methyl]-6-oxo-1H-purin-2-yl]acetamide; 9-[(2-Hydroxyethoxy)methyl]-N2-acetylguanineMolecular Formula: C₁₀H₁₃N₅O₄

Molecular weight: 267.24 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off-White Solid	-----
2	Purity by HPLC	>95.00 %	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*