


CERTIFICATE OF ANALYSIS

Reference Substance

Date

Simvastatin Lovastatin(Impurity E)
Batch No.

IUPAC Name: (1S,3R,7S,8S,8aR)-8-{2-[(2R,4R)-4-hydroxy-6-oxooxan-2-yl]ethyl}-3,7-dimethyl-1,2,3,7,8,8a-hexahydronaphthalen-1-yl (2S)-2-methylbutanoate


 Molecular Formula: C₂₄H₃₆O₅

Molecular weight: 404.54

Sr.#	Test	Results	Data ref.
1	Appearance	white colored powder	-----
2	Purity by HPLC	99.92%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C and protect from direct sunlight and moisture.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Date of manufacturing:
Retest date:
Document data reference:
Prepared by
Reviewed by
Approved by