

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Hydroxy Morpholino Thiadiazole

Batch No.

IUPAC Name:3-hydroxy 4 morpholino-1,2,5 thiadiazole

Molecular Formula:C₆H₉N₃O₂S

Molecular weight: 187.22 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	99.21%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at Refrigerator, Protect from direct sunlight**Shipment date:**

This certificate of analysis is valid for one year from date of shipment, provided the substance is stored under suitable conditions.

Document data reference:*Prepared by**Reviewed by**Approved by*