

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Tibolone

Batch No.

IUPAC Name: 17-Hydroxy-7a-methyl-19-nor-17a-pregn-5(10)-en-20-yn-3-one

Molecular Formula: C₂₁H₂₈O₂

Molecular weight: 312.45 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	97.21%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C and Protect from direct sunlight and moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of Manufacturing, provided the substance is stored under suitable conditions.

Manufacturing Date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*