


## CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

**Valaciclovir Related Compound G****Catalogue No.: VL149005****Batch No.:****IUPAC Name:** N,N-Dimethylpyridin-4-amineMolecular Formula: C<sub>7</sub>H<sub>10</sub>N<sub>2</sub>

Molecular Weight: 122.17 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	99.99%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

**Storage condition:** Store at 2-8 degree, and protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

**Manufacturing date:****Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*