

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date:

N-Formyl Valaciclovir

Catalogue no.: VL149006

Batch No.:

IUPAC Name: N-Formyl-L-valine 2-[(2-amino-1,6-dihydro-6-oxo-9H-purin-9-yl)methoxy]ethyl ester



Molecular Formula: C₁₄H₂₀N₆O₅

Molecular Weight: 352.35 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	96.71%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, and protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by