

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date:

9-Hydroxy Risperidone/Risperidone EP Impurity C**Catalog NO.: VL1520003**

Batch No.:

IUPAC Name: (9RS)-3-[2-[4-(6-Fluoro-1,2-benzisoxazol-3-yl)piperidin-1-yl]ethyl]-9-hydroxy-2-methyl-6,7,8,9-tetrahydro-4H-pyrido[1,2-a]pyrimidin-4-one

Molecular Formula: C₂₃H₂₇FN₄O₃

Molecular Weight: 426.48 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	white solid	-----
2	Purity by HPLC	98.80%	Chromatogram attached
3	MASS	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight & moisture.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:**

Document data reference:

*Prepared by**Reviewed by**Approved by*