

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Valsartan methyl Ester

Batch No.:

IUPAC Name: N-(1-Oxopentyl)-N-[[2'-(1H-tetrazol-5-yl)[1,1'-biphenyl]-4-yl]methyl]-L-valine methyl ester


 Molecular Formula: C₂₅H₃₁N₅O₃

Molecular Weight: 449.55 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White crystalline solid	-----
2	Purity by HPLC	99.51%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight and moisture.

Shipment date:

This certificate of analysis is valid for two years from date of Manufacturing, provided the substance is stored under suitable conditions.

Manufacturing Date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by