

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date:

Alosetron HCL Impurity D**Catalog NO: VL1730001**

Batch No.:

IUPAC Name: 2,5-dimethyl-2,3,4,5-tetrahydro-1H-pyrido[4,3-b]indol-1-one



Molecular Weight: 214.26g/mol

Molecular Formula: C₁₃H₁₄N₂O

Sr.#	Test	Results	Data ref.
1	Appearance	Cream colour solid	-----
2	Purity by HPLC	>95.00%	Chromatogram attached
3	Mass	Conform	Report attached
4	1H NMR	Conform	Report attached

Material Shipping Condition: At normal temperature**Long term Storage:** Store at 2-8°C, protect from moisture and direct sunlight.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*