


CERTIFICATE OF ANALYSIS

Reference Substance
Olmesartan Trityl Alcohol
Batch No.:
IUPAC Name: triphenylmethanol

Date



Molecular Weight: 260.33
Molecular Formula: C₁₉H₁₆O

Sr.#	Test	Results	Data ref.
1	Appearance	White solid Powder	-----
2	Purity by HPLC	99.9%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at refrigerator, protect from direct sunlight.

For Long term storage: 2-8°C,

Shipment date:

This certificate of analysis is valid for two year from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by