


## CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date

**Olmesartan EP impurity D/ Trityl Olmesartan****Catalogue no.: VL184005****Batch No.:****IUPAC Name:** ((5-methyl-2-oxo-1,3-dioxol-4-yl)methyl)4-(2-hydroxypropan-2-yl)-2-propyl-1((2'-(1-trityl-1H-tetrazol-5-yl)biphenyl-4-yl)methyl-1H-imidazole-5-carboxylate

Molecular Weight: 800.9

Molecular Formula: C<sub>48</sub>H<sub>44</sub>N<sub>6</sub>O<sub>6</sub>

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid powder	-----
2	Purity by HPLC	94.0%	Chromatogram attached
3	NMR	Conform	Report attached
4	Mass	Conform	Report attached

**Storage condition:** Store at refrigerator, protect from direct sunlight.

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

**Date of manufacturing:****Retest Date:**

Document data reference:

*Prepared by**Reviewed by**Approved by*