


## CERTIFICATE OF ANALYSIS

Reference Substance

Date

**Olmesartan Imp C (Olmesartan olefinic impurity)**
**Batch No.:**
**IUPAC Name:** (5-methyl-2-oxo-2H-1,3-dioxol-4-yl)methyl 4-(2-propen-2-yl)-2-propyl-1-({4-[2-(2H-1,2,3,4-tetrazol-5-yl)phenyl]phenyl}methyl)-1H-imidazole-5-carboxylate


Molecular Weight: 540.6

 Molecular Formula: C<sub>29</sub>H<sub>28</sub>N<sub>6</sub>O<sub>5</sub>

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid powder	-----
2	Purity by HPLC	98.5%	Chromatogram attached
3	NMR	Conform	Report attached
4	Mass	Conform	Report attached

**Storage condition:** Store at refrigerator, protect from direct sunlight.

**For Long term storage:** 2-8°C

**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

**Date of manufacturing:**
**Retest Date:**
**Document data reference:**
*Prepared by*
*Reviewed by*
*Approved by*