

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Azaerythromycin A (Azithromycin Impurity A)

Batch No.:

IUPAC Name: 6-demethyl Azithromycin



Molecular Formula: C₃₇H₄₀N₂O₁₂

Molecular Weight: 734.95g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White Solid	-----
2	Purity by HPLC	98.83%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by