

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Azithromycin Impurity C (3''-O-demethylazithromycin)

Batch No.

IUPAC Name: 3'-N-{[4-(Acetylamino)phenyl]sulfonyl}-3',3'-didemethylazithromycin



Molecular Weight: 734.98

Molecular Formula: C₃₇H₇₀N₂O₁₂

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	>95.0%	Chromatogram attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

Shipment Date:

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by