



## CERTIFICATE OF ANALYSIS

Reference Substance

Date

**Didecylamine HCl (Colsevelam impurity- 6)**

**Batch No. IUPAC Name:** Di- decylamine hydrochloride, N-decyldecane-1-amine hydrochloride



Molecular Weight: 298.2(Free Base)

:333.5(Salt)

Molecular Formula: C<sub>20</sub>H<sub>44</sub>ClN

Sr.#	Test	Results	Data ref.
1	Appearance	white solid	-----
2	Purity by IC	100.0%	-----
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

**Storage condition:** Store at refrigerator, protect from direct sunlight.

**For Long term storage:** 2-8°C,

**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

**Date of manufacturing:**

**Retest Date:**

**Document data reference:**

*Prepared by*

*Reviewed by*

*Approved by*