

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date:

Amlodipine Impurity 18**Catalog NO: VL19009**

Batch No.:

IUPAC Name: 2-((2-aminoethoxy)methyl)-4-(2-chlorophenyl)-1,4-dihydro-6-methylpyridine-3,5-dicarboxylic acid



Molecular Weight: 366.80 g/mol

Molecular Formula: C₁₇H₁₉ClN₂O₅

Sr.#	Test	Results	Data ref.
1	Appearance	Brown color liquid	-----
2	Purity by HPLC	86.17%	Chromatogram attached
3	Mass	Conform	Report Attached
4	¹ H NMR	Conform	Report Attached

Storage condition: Store at 2-8°C, protect from moisture and direct sunlight.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*