


CERTIFICATE OF ANALYSIS

Reference Substance

Date

Sertraline Impurity -F**Batch No.:**

IUPAC Name: (R)-4-(3,4-Dichlorophenyl)-1-tetralone

 $C_{16}H_{12}Cl_2O$
Mol. Wt.: 291.17

Molecular Weight: 291.17

Molecular Formula: $C_{16}H_{12}O$

Sr.#	Test	Results	Data ref.
1	Appearance	White colored powder	-----
2	Purity by HPLC	99.96%	Chromatogram attached
3	NMR	Conform	Report attached
4	Mass	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*