


CERTIFICATE OF ANALYSIS

Reference Substance

Date

Sertraline Related Compound-2**Batch No.**

IUPAC Name: (1R,4R)-Trans-4-(3,4-dichlorophenyl)-1,2,3,4-tetrahydro-N-methyl-1-naphthalenamine hydrochloride

Molecular Formula: C₁₇H₁₇Cl₂N

Molecular weight: 306.23

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	99.98%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight and moisture.**Shipment date:**

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*