

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date:

Sitagliptin Impurity-B

Catalogue no.: VL201006

Batch No.:

IUPAC Name:- (R)-3-amino-1-(3-(trifluoromethyl)-5,6-dihydro-[1,2,4]triazolo[4,3-a]pyrazin-7(8H)-yl)-4-(2,5-difluorophenyl)butan-1-one



Molecular Weight: 389.32g/mol

Molecular Formula: C₁₆H₁₆F₅N₅O

Sr.#	Test	Results	Data ref.
1	Appearance	Brown Sticky Mass	-----
2	Purity by HPLC	>90%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from moisture and direct sunlight.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by