

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date:

Mirtazapine N-Oxide /Mitrazapine EP Impurity-A/Mitrazapine Related Compound B

Catalogue no.: VL206004

Batch No.:

IUPAC Name: 1,2,3,4,10,14b-Hexahydro-2-methylpyrazino[2,1-a]pyrido[2,3-c][2]benzazepine 2-oxide monohydrate



Molecular Formula: C₁₇H₁₉N₃O

Molecular Weight: 281.35g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white Solid	-----
2	Purity by HPLC	95.79%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached
5	IR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by