

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard


Date:

1-keto Mirtazapine Impurity

Catalogue no.: VL206005

Batch No.:

IUPAC Name: 2-Methyl-3,4,10,14b-tetrahydrobenzo[c]pyrazino[1,2-a]pyrido[3,2-f]azepin-



Molecular Formula: C₁₇H₁₇N₃O

Molecular Weight: 279.3g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white Solid	-----
2	Purity by HPLC	98.56%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached
5	IR	Conform	Report attached

Storage condition: Store at 2-8 degree, protect from direct sunlight and moisture.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by