


CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD

Reference Standard

Date:

Losartan Impurity B**Catalogue No.:****Batch No.:****IUPAC Name:** (2'-(1H-tetrazol-5-yl)-[1,1'-biphenyl]-4-yl)methanol.Molecular Formula: C₁₄H₁₂N₄O

Molecular Weight: 252.27

Sr.#	Test	Results	Data ref.
1	Appearance	Light brown liquid	-----
2	Purity by HPLC	>95%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree and protect from moisture and Direct sunlight.

This certificate of analysis is valid for three years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*