

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date:

Losartan Bromo Impurity**Catalog NO: VL214010**

Batch No.:

IUPAC Name: 5-(4'-(bromomethyl)-[1,1'-biphenyl]-2-yl)-1H-tetrazole.

Molecular Formula: C₁₄H₁₁BrN₄

Molecular Weight: 315.17 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Pale yellow solid	-----
2	Purity by HPLC	89.18%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Note: Material is highly hygroscopic; open under flow of nitrogen.**Storage condition:** Store at 2-8°C and protect from moisture and direct sunlight.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*