

CERTIFICATE OF ANALYSIS


Reference Substance

Date:

Oxybutynin Impurity A (Cyclohexenyl analog of Oxybutynin Chloride)

Batch No.:

IUPAC Name: 4-(diethylamino)but-2-ynyl (RS)-2-(cyclohex-3-enyl)-2-cyclohexyl-2-hydroxyacetate.



Molecular Weight: 361.52

Molecular Formula: C₂₂H₃₅NO₃

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	91.48%	Chromatogram attached
3	Mass	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from moisture and direct sunlight.

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by