

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Ethylpropyl Analog/ (Oxybutynin Impurity E)

Batch No.:

IUPAC Name: 4-(ethylpropylamino)but-2-ynyl (±)-2-cyclohexyl-2-hydroxy-2-phenylacetate



Molecular Formula: C₂₃H₃₃NO₃

Molecular weight: 371.51 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	White solid	-----
2	Purity by HPLC	99.30%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at Refrigerator, Protect from direct sunlight.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest date:

Document data reference:

Prepared by

Reviewed by

Approved by