


CERTIFICATE OF ANALYSIS

Reference Substance

Date

3-Chloro Glyburide Impurity/ (Glibenclamide Impurity E)
Batch No.:
IUPAC Name: 1-cyclohexyl-3-[[4- [2- (3,5-dichloro-2-methoxybenzoyl)aminol-ethylphenyl]sulfonyl]urea

[3,5-dichloro-N-(4-(N-(cyclohexylcarbamoyl)sulfamoyl)phenethyl)-2-methoxybenzamide]


Molecular Formula: C₂₃H₂₇Cl₂N₃O₅S

Molecular weight: 528 g/mol

[M]⁻ m/z=526.6

Sr.#	Test	Results	Data ref.
1	Appearance	White to Off white solid	-----
2	Purity by HPLC	93.91%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform[M] ⁻ m/z=526.6	Report attached

Storage condition: Store at Refrigerator, Protect from moisture.

Shipment date:

This certificate of analysis is valid for one year from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:
Retest date:
Document data reference:

Prepared by

Reviewed by

Approved by