

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Paliperidone Impurity

Batch No.

IUPAC Name: 3-ethyl)-6,7,8,9-tetrahydro-9-hydroxy-2-methyl pyrido[1,2-a] pyrimidin-4-one

Molecular Formula: C₁₁H₁₆N₂O₂

Molecular weight: 208.26

Sr.#	Test	Results	Data ref.
1	Appearance	Pale yellow solid	-----
2	Purity by HPLC	97.01%	Chromatogram attached
3	NMR	Conform	Report attached
4	MASS	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest date:****Document data reference:***Prepared by**Reviewed by**Approved by*