


CERTIFICATE OF ANALYSIS

Reference Substance

Date

Rizatriptan impurity- C**Batch No.:****IUPAC Name:** N,N-Dimethyl-5-(1H-1,2,4-triazol-1-ylmethyl)-1H-indole-2-ethanamineMolecular Formula: C₁₅H₁₉N₅

Molecular Weight: 269.34

Sr.#	Test	Results	Data ref.
1	Appearance	Redish Brown residue	-----
2	Purity by HPLC	96.67%	Chromatogram attached
3	MASS	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight and moisture.**Shipment date:**

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*