

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Rizatriptan Desmethyl Impurity

Batch No.

IUPAC Name: 2-{5-[(1H-1,2,4-Triazol-1-yl)methyl]-1H-indol-3-yl}-N-methylethanamine



Molecular Weight: 255.32g/mol

Molecular Formula: C₁₄H₁₇N₅

Sr.#	Test	Results	Data ref.
1	Appearance	Dark brown solid	-----
2	Purity By HPLC	97.18%	Report attached
3	Mass	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*