


CERTIFICATE OF ANALYSIS

Reference Substance
Trospium Chloride Impurity D
 Batch No.:

Date:



Molecular Formula: C₂₈H₂₀O₄
 Molecular Weight: 420.46 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white colored powder	-----
2	Purity by HPLC	96.41%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree. Protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest Date:

Document data reference:

Prepared by

Reviewed by

Approved by