

CERTIFICATE OF ANALYSIS


Reference Substance

Date

Ziprasidone Open ring Impurity

Batch No.

IUPAC Name: [2-Amino-5-[2-[4-(1,2-benzisothiazol-3-yl)piperazin-1-yl]ethyl]-4-chlorophenyl]acetic Acid Sodium Salt


 Molecular Formula: C₂₁H₂₂ClN₄NaO₂S

Molecular weight: 452.93g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	97.69%	Chromatogram attached
3	MASS	Conform	Report attached
4	¹ H NMR	Conform	Report attached

Storage condition: Store at 2-8°C, protect from direct sunlight and moisture.

Shipment date:

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:

Retest date:

Document data reference:

Prepared by

Reviewed by

Approved by