

CERTIFICATE OF ANALYSIS OF REFERENCE STANDARD


Reference Standard

Date

Buspirone Related Compound N**Catalogue No.: VL249006**

Batch No.:

IUPAC Name: 8, 8'-(butane-1, 4-diyl) bis (8-azaspiro [4, 5] decane-7, 9-dione).

Molecular Formula: C₂₂H₃₂N₂O₄

Molecular Weight: 388.50 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Off white solid	-----
2	Purity by HPLC	91.16%	Chromatogram attached
3	MASS	Conform	Report attached
4	NMR	Conform	Report attached

Storage condition: Store at 2-8°C, Protect from direct sunlight & moisture.

This certificate of analysis is valid for three years from date of manufacturing provided the substance is stored under suitable conditions.

Date of manufacturing:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*