

CERTIFICATE OF ANALYSIS

Reference Substance


Date:

4-chlorobupropion**Batch No.:**

IUPAC Name: 2-(tert-Butylamino)-1-(3,4-dichlorophenyl)propan-1-one

OR

2-(tert-butylamino)-3'-4'-chloropropiophenone hydrochloride

Molecular Formula: C₁₃H₁₇Cl₂NO

Molecular Weight: 274.0 g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	A Cream colour solid	-----
2	Purity by HPLC	99.69%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, and protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*