

CERTIFICATE OF ANALYSIS

Reference Substance

Date:


Bupropion Related Compound C

Batch No.:

IUPAC Name: 1-(3-Chlorophenyl)-2-hydroxy-1-propanone

OR

1-(3-Chlorophenyl)-2-hydroxypropan-1-one

Molecular Formula: C₉H₉ClO₂

Molecular Weight: 184.62g/mol

Sr.#	Test	Results	Data ref.
1	Appearance	Pale yellow liquid	-----
2	Purity by HPLC	95.08%	Chromatogram attached
3	Mass	Conform	Report attached
4	H NMR	Conform	Report attached

Storage condition: Store at 2-8 degree, and protect from direct sunlight and moisture.

This certificate of analysis is valid for two years from date of manufacturing, provided the substance is stored under suitable conditions.

Manufacturing date:**Retest Date:****Document data reference:***Prepared by**Reviewed by**Approved by*